

1.67 Clear & Photochromic FT-28, FT-35 & 7x28 Semi-Finished

clear


	NEW		
FT-28 SF Clear	FT-35 SF Clear		7x28 SF Clear
Base: 1, 3, 5, 7	Base: .1, 3, 5, 7		Base: 1, 3, 5, 7
Adds: 1.00D to 3.00D	Adds: 1.00D to 3.50D		Adds: 1.50D to 3.50D
Diameter: 75mm	Diameter: 70mm		Diameter: 70mm
Material: MR-7	Material: MR-10		Material: MR-10

sunSMART^{II}
grey™

sunSMART^{II}
brown™

ultraSUN^{II}
sunglass photo

FT-28 SF Photo
Base: .50, 2, 4, 6, 8
Adds: 1.00D to 3.00D
Diameter: 73mm
Material: MR-10



7x28 SF Photo
Base: 2, 4, 6
Adds: 1.50D to 3.50D
Diameter: 75mm
Material: MR-10



sunSMART^{II}
grey™

- Premium photochromics and clear to complete your 1.67 multifocal product line.
- sunSMART photochromics with enhanced color performance, high-temperature stability, and enduring color consistency.
- ultraSUN sunglass photochromic provides a sunglass dark photochromic lens option when activated.
-  Photochromics developed and blended in the USA.