



1.67 Clear & Photochromic

FT-28 & 7x28 Semi-Finished

Premium photochromics and clear to complete your 1.67 multifocal product line.

Photochromics with enhanced color performance, high-temperature stability, and enduring color consistency.



Photochromics developed and blended in the USA.

sunSMART II
grey™

Availability
Base: .50, 2, 4, 6, 8
Adds: 1.00D to 3.00D
Diameter: 73mm
Material: MR-10

sunSMART II
brown™

ultraSUN II
sunglass photo™




FT-28 Semi- Finished

clear

Base: 1, 3, 5, 7
Adds: 1.00D to 3.00D
Diameter: 75mm
Material: MR-7

sunSMART II
grey™

Availability
Base: 2, 4, 6
Adds: 1.50D to 3.50D
Diameter: 75mm
Material: MR-10



7x28 Semi- Finished

clear

Base: 1, 3, 5, 7
Adds: 1.50D to 3.50D
Diameter: 70mm
Material: MR-10



ultraSUN II
sunglass photo™
Now Available!
Sunglass dark when activated with exceptional high temperature performance. Perfect for an active outdoor lifestyle.