

# 1.67 & 1.74 Clear & Photochromic FT-28, FT-35 & 7x28 Semi-Finished

clear

<b>FT-28 SF Clear</b>
Base: 1, 3, 5, 7
Adds: 1.00D to 3.00D
Diameter: 75mm
Material: MR-7

<b>FT-35 SF Clear</b>
Base: .1, 3, 4, 7
Adds: 1.50D to 3.50D
Diameter: 70mm
Material: MR-10


<b>7x28 SF Clear</b>
Base: 1, 3, 5, 7
Adds: 1.50D to 3.50D
Diameter: 70mm
Material: MR-10

sunSMART II grey™

sunSMART II brown™

ultrasun II sunglass photo™


<b>FT-28 SF Photo</b>
Base: .50, 2, 4, 6, 8
Adds: 1.00D to 3.00D
Diameter: 73mm
Material: MR-10



sunSMART II grey™


**NEW LENS 1.74**

<b>FT-28 SF Photo</b>
Base: .50, 1.50, 3.50, 8.50
Adds: 1.00D to 3.00D
Diameter: 70,75,75,70mm
<b>Material: 1.74</b>



sunSMART II grey™

<b>7x28 SF Photo</b>
Base: 2, 4, 6
Adds: 1.50D to 3.50D
Diameter: 75mm
Material: MR-10

- Premium photochromics and clear to complete your 1.67 & 1.74 multifocal product line.
- sunSMART photochromics with enhanced color performance, high-temperature stability, and enduring color consistency.
- ultrasun sunglass photochromic provides a sunglass dark photochromic lens option when activated.
-  Photochromics developed and blended in the USA.